

244 Bariatric

Barrier Free Exam Table



Healthcare Furniture

244 Bariatric

Barrier Free Exam Table



The Ritter 244 Bariatric Power Treatment Table is designed to provide accessibility and dignity for patients while supporting safe, efficient care for clinical staff. With an 850 lb weight capacity, the table is built to safely lift special-needs patients repeatedly and reliably.

Depth	80"
Width	32"
Height	18-34"
Seat Width	32"
Seat Height	18-34"
Weight Capacity	850 lbs

244 Bariatric

Barrier Free Exam Table



OPTIONS

