

500 series

Hydraulic Stools



Healthcare Furniture

500 series

Hydraulic Stools



The 500 Series Hydraulic Stool features a solid five-leg base that extends beyond the seat diameter for maximum stability in any position. Its 2.5-inch casters are available in either hard or soft rubber to suit the needs of the space. The stool is offered in standard, medium, or high-height variants, along with four seat styles.



540

Round Stool

Seat Size	15"
Base Diameter	19.75"
Low Height	17.5-25"
Standard Height	20.5-28"
High Height	24-31.5"
Weight Capacity	300 lbs



546

Round Stool with Back

Seat Size	15"
Base Diameter	19.75"
Low Height	17.5-25"
Standard Height	20.5-28"
High Height	24-31.5"
Weight Capacity	300 lbs



548

Rectangle Seat with Back

Seat Size	14 x 16"
Base Diameter	19.75"
Low Height	17.5-25"
Standard Height	20.5-28"
High Height	24-31.5"
Weight Capacity	300 lbs



556

Round Stool with Arm

Seat Size	15"
Base Diameter	19.75"
Low Height	17.5-25"
Standard Height	20.5-28"
High Height	24-31.5"
Weight Capacity	300 lbs