

600 Series Operators Stool



Healthcare Furniture

600 Series Operators Stool



Designed to support a wide range of sitting positions, the uniquely shaped backrest allows for greater freedom of movement, variation, and natural posture. The saddle seat adjusts from low to high working positions, including standing height, helping users move more dynamically throughout the day. Often referred to as the “Sonographer’s Chair,” its backrest provides full trunk and forearm support when seated backwards and arm support when seated sideways.

Depth	28.75"
Width	28.75"
Height	34.5-39.5"
Seat Depth	16.5"
Seat Width	18.5"
Standard Seat Height	16.5-21.5"
Medium Seat Height	19-26"
High Seat Height	22.5-32"

600 Series Operators Stool



OPTIONS

