

Bariatric Phlebotomy Chair



Healthcare Furniture

Bariatric Phlebotomy Chair



The PC8690 Bariatric Blood Draw Chair features seamless upholstery and a heavy-duty steel frame, providing a 750 lb weight capacity. Its 26.625-inch-wide seat accommodates both bariatric patients and parents with small children, while the flip-up, height-adjustable arm board helps position the arm for optimal support.

Depth	35"
Width	41"
Height	41.75"
Seat Width	26.625"
Seat Height	22"
Weight Capacity	750 lbs