

Cama

Healthcare Recliner

 **statmedical**
we take care



Healthcare Furniture

Cama

Healthcare Recliner



The Cama Care Recliner combines smooth reclining motion with the therapeutic styling of the Cama chair series. Independent foot and backrest controls allow for flexible positioning from fully upright to a flat sleeping position. A composite foam seat with dynamic suspension provides lasting comfort, while the crumb-space clean-out makes cleaning and maintenance easier. Comes standard with polyester arm caps for protection and longevity.



108-56

Infinite Locking Recliner

Depth	37"
Width	32"
Height	45"
Seat Depth	20"
Seat Width	24"
Seat Height	18"
Arm Height	25"
Reclined Depth	72"
Weight Capacity	550 lbs

Cama

Healthcare Recliner



OPTIONS



Drainage Bag Hook-up



Oxygen Tank Holder



Seat Moisture Barrier



Footrest guards



Folding Arm



Rolling Foot Tray



IV Pole and Holder



Heated Seat



Infection Control Flip-Up Seat

Cama

Healthcare Recliner



Soft Grip Push Bar



Solid Surface Arm Caps



Swivel Profile Tablet



Folding Side Tablet



Swiveling Side Tablet



Power/USB



Central Lock



Restraint Hooks