

# Carson Recliner



# Comfort and Functionality in One Piece

The Carson Recliner is a contemporary solution for demanding healthcare environments, including triage areas, infusion treatment spaces, and patient rooms. With a clean, modern aesthetic, it is designed to provide a high level of comfort for patients while ensuring ease of use for caregivers and family members.



## Features

- Available in 23" and 26" seat widths
- Field-replaceable upholstery covers
- Back construction features 14-gauge inner steel structure
- Lever positioning permits both caregiver and patient to operate footrest
- Ergonomically-positioned handles for easy transport
- Side-folding transfer arm

## Options:

- Central locking casters for mobility and stability
- Side tablets for convenient access to essentials
- Slide-out foot tray for added functionality
- Fold-down arm for versatile usage
- Trendelenburg's position for enhanced patient care
- IV pole attachment for medical procedures
- Foley bag hook for patient convenience
- Moisture barrier for durability

	Depth	Width	Seat Width	Height	Seat Height	Weight capacity
CRR123F	36"	29"	23"	48"	20.5"	500 lbs
CRR126F	36"	32"	26"	48"	20.5"	500 lbs
CRR123FT	36"	29-32"	23"	48"	20.5"	500 lbs
CRR126FT	36"	32-35"	26"	48"	20.5"	500 lbs



CRR123BH



CRR126F



CRR123FT



CRR126FT