

Carson Recliner



Healthcare Furniture

Carson

Recliner



The Carson Recliner combines a clean, modern aesthetic with thoughtful design for demanding healthcare environments, including triage areas, infusion spaces, and patient rooms. It provides a high level of patient comfort while remaining easy for caregivers and family members to use.



CRR123F
23" Recliner



CRR126F
26" Recliner



CRR123FT
23" Recliner, Patient Transfer Arm



CRR126FT
26" Recliner, Patient Transfer Arm

CRR123F 23" Recliner		CRR126F 26" Recliner		CRR123FT 23" Recliner, Patient Transfer Arm		CRR126FT 26" Recliner, Patient Transfer Arm	
Depth	36"	Depth	36"	Depth	36"	Depth	36"
Width	29"	Width	32"	Width	29/32"	Width	32/35"
Height	48"	Height	48"	Height	48"	Height	48"
Seat Depth	20"	Seat Depth	20"	Seat Depth	20"	Seat Depth	20"
Seat Width	23"	Seat Width	26"	Seat Width	23"	Seat Width	26"
Seat Height	20.5"	Seat Height	20.5"	Seat Height	20.5"	Seat Height	20.5"
Arm Height	26"	Arm Height	26"	Arm Height	26"	Arm Height	26"
Reclined Depth	75"	Reclined Depth	75"	Reclined Depth	75"	Reclined Depth	75"
Weight Capacity	500 lbs	Weight Capacity	500 lbs	Weight Capacity	500 lbs	Weight Capacity	500 lbs

Carson

Recliner

OPTIONS



Moisture Barrier



Pedal Central-Locking Casters



IV Pole



Trendelenburg Position



Accessory Hook



Central Locking Casters



Flip-Swivel Table



Flip-Up Side Table



Polyurethane Arm Caps

Carson Recliner

OPTIONS



Wood Arm Caps



Solid Surface Arm Caps

Solid Surface



Glacier White



Linen



Pearl Grey



Silver Grey



Bisque



Cameo White



City Roast



Dove

Carson

Recliner

Wood



Aria



Candlelight



Chestnut



Chocolate



Dolce Vita



Grenada



Hardrock Maple



Havana



Palomino



Pebble




Sand



White Chocolate

Polyurethane



Black



Taupe



Grey

