

Chrome

Step Stool



Minor Equipment

Chrome Step Stool



This chrome stool features a fully welded steel base and top with no screws or rivets for added strength and durability. A rubber tread surface provides secure footing, while reinforced rubber feet offer stability. Finished in durable chrome, the stool is available with or without a handrail to suit different needs.



SSC-221

Chrome Step Stool

Depth	11.375"
Width	14.25"
Platform Height	9"
Weight Capacity	350 lbs



SSC-221H

Chrome Step Stool with Handrail

Depth	11.375"
Width	15.25"
Platform Height	9"
Handle Height	44.75"
Weight Capacity	350 lbs