

Clifton

Stainless Steel Stool



Healthcare Furniture

Clifton

Stainless Steel Stool



The Clifton Stainless Steel Stool features a recessed 15-inch seat that adjusts easily by turning and includes a safety stop to prevent inadvertent removal. The polished seat is designed to accept a vinyl insert or seat pad, while a circular heel rest and flared legs add practical support and stability.



7745SS

Stool, Rubber Feet

Seat Diameter	15"
Height	19-31"
Weight Capacity	300 lbs



7745SS-C

Stool, Casters

Seat Diameter	15"
Height	21-33"
Weight Capacity	300 lbs



7745SS-P

Upholstered Stool, Rubber Feet

Seat Diameter	15"
Height	19-31"
Weight Capacity	300 lbs



7745SS-PC

Upholstered Stool, Casters

Seat Diameter	15"
Height	21-33"
Weight Capacity	300 lbs