

Comfort-3 Eco

Powered Patient Treatment Chair



Healthcare Furniture

Comfort-3 Eco

Powered Patient Treatment Chair



The Comfort-3 Eco Chair combines modern design, exceptional patient comfort, and practical space-saving efficiency. Ideal for dialysis, infusion therapy, chemotherapy, and blood collection applications, its compact footprint can free up to 15% more clinical floor space while providing easy front entry and exit for patients. Standard powered positioning makes it simple for either staff or patients to change the chair from upright to horizontally flat, or into trendelenburg position.



Comfort-3 Eco
Powered Patient Treatment Chair

Depth 39"
Width 36.5"
Height 54.75"

Seat Width 23.75"
Seat Height 21.75"
Weight Capacity 440 lbs

Backrest Adjustment 80° to -14°
Seat Tilt Adjustment -14°
Legrest Adjustment 0° to 90°

Comfort-3 Eco

Powered Patient Treatment Chair



OPTIONS



EPN Kit




Trendelenburg Pedals



Battery Backup



Power/USB



Multifunction Table



Adjustable Arms



Folding Side Rails




Plastic Side Rails



Push Bar



Push Bar Pro



Phone/Tablet Clamp



Seat Side Handles

Comfort-3 Eco

Powered Patient Treatment Chair

OPTIONS

