

Connect

Electric Recliner



Healthcare Furniture

Connect

Electric Recliner



The Connect Electric Exam Recliner combines reclining, exam table functionality, and safe patient handling for lateral transfers. Its height-adjustable design raises to a comfortable working height for caregivers, while optional accessories expand its use as an exam table. A continuous footrest comes standard for easier cleaning and reduced wear on the footplates, making Connect a versatile solution for acute care environments and beyond.



Connect-E29
Recliner

Depth	36"
Width	30.5"
Height	47.5"
Seat Depth	22"
Seat Width	24.5"
Seat Height	21-32"
Reclined Depth	73.5"
Arm Height	24"
Weight Capacity	650 lbs



Connect-E29-T
Recliner, Trendelenburg

Depth	38"
Width	30.5"
Height	47.5"
Seat Depth	22"
Seat Width	24.5"
Seat Height	21-32"
Reclined Depth	73.5"
Arm Height	24"
Weight Capacity	650 lbs

Connect

Electric Recliner

OPTIONS



Moisture Barrier



Pedal Central-Locking Casters



Toggle Central-Locking Casters



Button Folding Arm



Drainage Bag Hook



Rolling Foot Tray



IV Pole



Footrest Guard



Fall Prevention Footrest

Connect

Electric Recliner

OPTIONS



Push Bar



Folding Side Table



Flip-Swivel Side Table



Back Style




Speaker Headrest



Heated Seat and Back



Massage Back



Ergotron Arm



Stirrups

Connect

Electric Recliner

Polyester



Wood

