

Connect

Hydraulic Recliner



Customizable by Design

The Connect Hydraulic Recliner was built with patient and caregiver comfort in mind. This endlessly versatile recliner can transform into an exam table with enough height to accommodate the caregiver. Its hydraulic mechanisms and sturdy design ensure safe transportation and adjustments. Connect offers many customization options that increase its applications.



Features

- 12-year warranty
- Hydraulic height adjustment
- Mobile transport base
- Lays completely flat
- Easy lateral transfers

Options

- Stirrups
- Folding or swivel side tablet
- Fold-down arm
- Rolling foot tray
- Soft grip push bar
- Ergotron arm
- Drainage bag holder
- Heated seat and massage
- Central locking casters
- Straight or flared back
- Trendelenburg
- Battery backup
- 3-button foot control
- IV pole
- Push bar
- Paper roll dispenser

	Depth	Width	Seat Width	Height	Seat Height	Weight Capacity
Connect-1	38"	34"	24"	47"	18" - 32"	500 lbs
Connect-2	36"	33"	24"	46"	18" - 32"	500 lbs



Connect-1



Connect-2