

Connect

Hydraulic Recliner



Healthcare Furniture

Connect

Hydraulic Recliner



The Connect combines the functions of a recliner, exam table, and safe patient handling solution for lateral transfers. Its height-adjustable design raises to a comfortable working height for caregivers, while optional accessories allow it to be customized for exam use. With broad applications across acute care facilities and beyond, Connect offers exceptional flexibility in one versatile solution.



Connect-1

Lift Recliner, Trendelenburg

Depth	38"
Width	34"
Height	47-61"
Seat Depth	22"
Seat Width	24"
Seat Height	18-32"
Reclined Depth	74"
Arm Height	27-41"
Weight Capacity	550 lbs



Connect-2

Lift Recliner

Depth	36"
Width	33"
Height	46-60"
Seat Depth	21"
Seat Width	24"
Seat Height	18-32"
Reclined Depth	74"
Arm Height	26-40"
Weight Capacity	550 lbs

Connect

Hydraulic Recliner

OPTIONS



Moisture Barrier



Pedal Central-Locking Casters



Toggle Central-Locking Casters



Button Folding Arm



Lever Folding Arm



Drainage Bag Hook



Rolling Foot Tray



IV Pole



Footrest Guard

Connect

Hydraulic Recliner

OPTIONS



Connect

Hydraulic Recliner

Polyester



Wood

