

Drop-Arm Patient Recliner



Healthcare Furniture

Drop-Arm Patient Recliner



The Drop Arm Recliner is designed to support patient transfers, comfort, and versatility in clinical environments. Its arms lower flush with the seat, while non-porous blow-molded fold-down trays are included on both sides. A standard Trendelenburg release is accessible from either side, and features such as an ergonomic headrest, adjustable elastic back supports, and an open elastic seat base enhance comfort and durability. Suitable for Dialysis, Oncology, Acute Care, and other clinical settings.



6550

Recliner

Depth	40"
Width	33.75"
Height	49"
Seat Depth	20"
Seat Width	21"
Seat Height	21"
Reclined Depth	76.5"
Arm Height	21-27"
Weight Capacity	275 lbs




6570

Bariatric Recliner


Depth	43"
Width	37.75"
Height	49"
Seat Depth	20"
Seat Width	25"
Seat Height	21"
Reclined Depth	76.5"
Arm Height	21-27"
Weight Capacity	450 lbs

Drop-Arm Patient Recliner

OPTIONS



Heated Seat



Heat and Massage Seat



Headrest Protector



IV Pole

