

Eclipse

Overbed Table



Healthcare Furniture

Eclipse

Overbed Table



The Eclipse Overbed Table is a budget-friendly option designed to provide dependable quality at an economical price. It features an attractive thermoformed top designed to contain large spills, along with a sturdy powder-coated H-base. A spring-assisted lift mechanism provides smooth, infinite height adjustment from 29 to 45 inches.

Depth	16"
Width	30"
Height	29-45"
Base Height	3.25"
Weight Capacity	25 lbs

Eclipse

Overbed Table

OPTIONS



Thermofoil




Aria



Candlelight



Chestnut



Chocolate



Dolce Vita



Grenada



Hardrock Maple



Havana



Palomino



Pebble



Sand



White Chocolate