

Guardian

Bassinets



Healthcare Furniture

Guardian Bassinet




The M-22 Bassinet combines durable construction with practical infant-care features. Trendelenburg and reverse-Trendelenburg positioning are available through the C-shaped hook system, while easy-to-clean stainless steel handles and KYDEX surfaces support durability and infection control. The clear PETG tub offers four times greater impact resistance than acrylic tubs, and the base panel helps reduce shin contact.

Depth	21.75"
Width	34"
Height	36"
Surface Height	30"

Guardian Bassinet

OPTIONS



Kydex



American
Natural



Aria



Blanca Elme



Bright White



Eiderdown



Summer Flame



Walnut