

H Series

Phlebotomy Chair



Healthcare Furniture

H Series

Phlebotomy Chair



The H Series blood drawing chairs feature a convenient hydraulic lift, adjusted using a large pump-action foot pedal to position patients at the ideal working height. Models 6310 and 6330-P include a single stationary arm and straight flip arm, while 6320 and 6340-P feature a drawer arm and straight flip arm. Choose between a comfortable vinyl-upholstered seat or the durable EZ-Clean single-piece solid plastic seat for improved cleanability.



6310

Padded Seat



6320

Padded Seat



6330-P

EZ-Clean Seat



6340-P

EZ-Clean Seat

Depth	30"	Depth	30"	Depth	30"	Depth	30"
Width	30.75"	Width	38"	Width	30.75"	Width	38"
Seat Width	21"	Seat Width	21"	Seat Width	18"	Seat Width	18"
Seat Height	20.25-27.5"	Seat Height	20.25-27.5"	Seat Height	20.25-27.5"	Seat Height	20.25-27.5"
Weight Capacity	300 lbs	Weight Capacity	300 lbs	Weight Capacity	300 lbs	Weight Capacity	300 lbs