

Ivy

Lounger Recliner

 **statmedical**
we take care



Healthcare Furniture

www.statmedicalcanada.com | 866 783 4076 | sales@statmedicalcanada.com

Ivy

Lounger Recliner



The Ivy Lounger provides a greater sense of personal space with its integrated wings, creating a comfortable and supportive seating experience. The chair extends into a relaxing lounge position for added comfort during longer stays, while its unique wall-saver design allows it to be positioned directly against the wall without interfering with the reclining function.



813-27W

Wingback Lounger

Depth	36"
Width	31"
Height	42"

Seat Depth	21"
Seat Width	24"
Seat Height	18"
Reclined Depth	72"
Arm Height	25"
Weight Capacity	350 lbs



813-27W-2T

Dual Wingback Lounger

Depth	36"
Width	77"
Height	42"

Seat Depth	21"
Seat Width	24"
Seat Height	18"
Reclined Depth	72"
Arm Height	25"
Weight Capacity	700 lbs

OPTIONS

