

# Kangaroo

## Healthcare Recliner



## Healthcare Furniture

# Kangaroo

## Healthcare Recliner



Kangaroo is a popular recliner in maternal newborn areas as a great option for nursing mothers, with increased lateral support and comfort from the padded back. These features also are a great option whenever departments are looking for a more comfortable option than traditional recliners. A variety of models allow for sleep functionality, soothing glide movement, and more.



**RC523-10**  
3 Position Recliner

|                |     |
|----------------|-----|
| Depth          | 39" |
| Reclined Depth | 70" |
| Width          | 29" |
| Height         | 44" |

|                 |         |
|-----------------|---------|
| Seat Depth      | 22"     |
| Seat Width      | 24"     |
| Seat Height     | 20"     |
| Weight Capacity | 350 lbs |



**RC523-15**  
Multi-Lock Recliner

|                |     |
|----------------|-----|
| Depth          | 36" |
| Reclined Depth | 68" |
| Width          | 33" |
| Height         | 48" |

|                 |         |
|-----------------|---------|
| Seat Depth      | 20"     |
| Seat Width      | 23"     |
| Seat Height     | 18"     |
| Weight Capacity | 350 lbs |



**RC523-35**  
Gliding Recliner

|                |     |
|----------------|-----|
| Depth          | 39" |
| Reclined Depth | 66" |
| Width          | 32" |
| Height         | 44" |

|                 |         |
|-----------------|---------|
| Seat Depth      | 20"     |
| Seat Width      | 24"     |
| Seat Height     | 19"     |
| Weight Capacity | 350 lbs |



**RC523-51**  
Multilock Sleeper

|                |     |
|----------------|-----|
| Depth          | 35" |
| Reclined Depth | 71" |
| Width          | 30" |
| Height         | 48" |

|                 |         |
|-----------------|---------|
| Seat Depth      | 20"     |
| Seat Width      | 22.5"   |
| Seat Height     | 19"     |
| Weight Capacity | 550 lbs |

# Kangaroo

## Healthcare Recliner



**RC523-52**

Sleeper, Trendelenburg

|                |     |
|----------------|-----|
| Depth          | 36" |
| Reclined Depth | 71" |
| Width          | 32" |
| Height         | 48" |

|                 |         |
|-----------------|---------|
| Seat Depth      | 19"     |
| Seat Width      | 23"     |
| Seat Height     | 18"     |
| Weight Capacity | 550 lbs |



**RC523-53**

Multilock Gliding Sleeper

|                |     |
|----------------|-----|
| Depth          | 36" |
| Reclined Depth | 71" |
| Width          | 32" |
| Height         | 48" |

|                 |         |
|-----------------|---------|
| Seat Depth      | 20"     |
| Seat Width      | 24"     |
| Seat Height     | 18"     |
| Weight Capacity | 550 lbs |




# Kangaroo

## Healthcare Recliner

### Polyester Arm Caps



### Thermofoil Tables



# Kangaroo

## Healthcare Recliner

### OPTIONS / ACCESORIES



# Kangaroo

## Healthcare Recliner

### OPTIONS / ACCESORIES



Folding Arm - Lever



Footrest Guards



Heated Seat



IPad Arm



IV Pole and Holder



Rolling Foot Tray



Soft Grip Push Bar