

L150

Lifedop doppler

 **statmedical**
we take care



Affordable, Reliable Quality

The LifeDop L150 is an accurate Doppler probe that has proven excellent for obstetrical and vascular practices. Its interchangeable probes make it a versatile addition to any practice, offering a host of capabilities. The LifeDop L150 operates off of 3 AA batteries that will last up to 500 assessments. The energy-conscious Doppler performs with a clarity of sound that compares to the most expensive models in the market.



Features

- AA Alkaline battery operation
- Made in America
- Energy saving mode
- Choose from a large family of probes

Options

- L150 – Standard
- L150A – Audio Recording
- L150R – Rechargeable

Probes

- SDW – 2MHz Waterproof Obstetrical Probe
- SD2 – 2MHz Late-Term Obstetrical Probe
- SD3 – 3MHz Early-Term Obstetrical Probe
- SD4 – 4MHz Deep Vascular
- SD5 – 5MHz Mid Vascular
- SD8 – 8MHz Peripheral Vascular

	Depth	Width	Height	Weight
Main Unit	2.75"	5.5"	1.4"	0.7 lbs

