

Lab Series

Phlebotomy Chair



Healthcare Furniture

Lab Series

Phlebotomy Chair



The Lab Series phlebotomy chairs feature a one-piece black plastic seat that is easy to clean, along with a flip arm that is adjustable in both height and depth. Designed for practical everyday use, the Lab Series is well suited for both large and small facilities.



66060-P

Phlebotomy Chair



66029-P

Phlebotomy Chair with Drawer



66011-P

Extra-Tall Phlebotomy Chair with
Drawer



66022-P

Extra-Tall Phlebotomy Chair

Depth	28"	Depth	28"	Depth	28"	Depth	28"
Width	31"	Width	40"	Width	31"	Width	40"
Seat Depth	17"	Seat Depth	17"	Seat Depth	17"	Seat Depth	17"
Seat Width	18"	Seat Width	18"	Seat Width	18"	Seat Width	18"
Seat Height	20"	Seat Height	20"	Seat Height	26"	Seat Height	26"
Arm Height	27-33"	Arm Height	27-33"	Arm Height	33-39"	Arm Height	33-39"
Weight Capacity	300 lbs	Weight Capacity	300 lbs	Weight Capacity	300 lbs	Weight Capacity	300 lbs