

Mayo Stand Chrome Economy



Minor Equipment

Mayo Stand

Chrome Economy



The SSC-126 Mayo Stand is an economical option designed for light-duty use. It features a removable stainless steel tray, a 1-inch chrome-plated pole, and a traditional sled-style base with two ball-bearing rubber casters. A screw-style thumb knob allows for easy height adjustment and secure positioning.

Depth	18.25"
Width	31.5"
Height	34.5-54"
Tray Depth	12.625"
Tray Width	19.125"
Load Capacity	15 lbs