

Mayo Stand Chrome Mobile



Minor Equipment

Mayo Stand Chrome Mobile



The SSC-124 Mayo Stand features a chrome-plated pole, stainless steel tray, and stable five-leg cast aluminum base with a wall-saver feature. Its fully mobile design allows for easy positioning, making it ideal for clinic and seated procedures. A traditional thumb knob provides secure height adjustment.

Depth	12.625"
Width	19.125"
Height	29.5-49"
Tray Depth	12.75"
Tray Width	19"
Weight Capacity	20 lbs