

Medlife

Exam Table



Healthcare Furniture

Medlife

Exam Table



The Medlife Exam Table is a budget-friendly option designed to meet the needs of medical facilities without compromising quality. It features seamless premium vinyl upholstery and a fully welded, powder-coated all-steel frame for durability and comfort. Available with either a front or side step stool, the table also includes pass-through rear drawers accessible from either side, an integral table paper holder, a pneumatically adjustable back, and a pull-out foot section that extends the table to 75 inches when laid flat.



2250X

Front Step, Black Upholstery



2250X

Front Step, Mocha Upholstery



2250X

Front Step, Blue Upholstery



3250X

Side Step, Black Upholstery



3250X

Side Step, Mocha Upholstery

Depth	58.5"	Depth	58.5"	Depth	58.5"	Depth	58.5"	Depth	58.5"
Width	27.5"	Width	27.5"	Width	27.5"	Width	27.5"	Width	27.5"
Height	34"	Height	34"	Height	34"	Height	34"	Height	34"
Extended Depth	75"	Extended Depth	75"	Extended Depth	75"	Extended Depth	75"	Extended Depth	75"
Step Stool Height	9.5"	Step Stool Height	9.5"	Step Stool Height	9.5"	Step Stool Height	9.5"	Step Stool Height	9.5"
Weight Capacity	500 lbs	Weight Capacity	500 lbs	Weight Capacity	500 lbs	Weight Capacity	500 lbs	Weight Capacity	500 lbs