

NC24C16-DT2-AVD

Narcotics Cabinet with Audio/Visual Alarm



Minor Equipment

NC24C16-DT2-AVD

Narcotics Cabinet with Audio/Visual Alarm



The NC24C16-DT2-AVD Narcotics Cabinet provides secure medication storage with dual locking doors and integrated visual and audible alarms. A flashing light and 84-decibel alarm activate when the cabinet is opened, providing added security. Its wall-mountable, vertical design features two adjustable or removable shelves for flexible storage of controlled substances.

Depth	8"
Width	16"
Height	26"