

Oasis

Treatment Recliner



Healthcare Furniture

Oasis

Treatment Recliner



The Oasis Treatment Recliner is thoughtfully designed for demanding 24/7 healthcare environments, combining durability with exceptional comfort.

Generous padding and a gas-cylinder mechanism allow for smooth, infinite recline positioning, helping create a more relaxing experience in clinical settings.



SR900-23

23" Seat

Depth	36"
Width	29"
Height	48"
Seat Depth	20"
Seat Width	23"
Seat Height	20.5"
Arm Height	27.25"
Reclined Depth	75"
Weight Capacity	500 lbs



SR900-26

26" Seat

Depth	36"
Width	32"
Height	48"
Seat Depth	20"
Seat Width	26"
Seat Height	20.5"
Arm Height	27.25"
Reclined Depth	75"
Weight Capacity	500 lbs

Oasis

Treatment Recliner

OPTIONS



Moisture Barrier




Foot Tray



IV Pole



Accessory Hook



Central Locking Casters



Flip-Swivel Table



Flip-Up Side Table



Trendelenburg Positioning



Polyurethane Caps

Oasis

Treatment Recliner

OPTIONS




Wood Arm Caps

Solid Surface




Glacier White



Linen




Pearl Grey



Silver Grey



Bisque



Cameo White



City Roast



Dove

Oasis

Treatment Recliner

Wood



Aria



Candlelight



Chestnut



Chocolate



Dolce Vita



Grenada



Hardrock Maple



Havana



Palomino



Pebble




Sand



White Chocolate

Polyurethane



Black



Taupe



Grey

