

P-Series

Stainless Steel Bassinet



Healthcare Furniture

P-Series

Stainless Steel Bassinet



The P-Series Bassinet is a premium stainless steel collection designed for practical infant care. The basket can be raised at either end, while the slim profile allows caregivers close access from any side. A foot-end drawer provides ample storage while remaining neatly tucked beneath the bassinet.



P-1110-SS

Basic Bassinet



P-1110-A-SS

Basic Bassinet, Front Drawer



P-1110-B-SS

Basic Bassinet, Side Drawer



P-1115-A-SS

Cabinet Base Bassinet



P-1115-SS

Cabinet Base Bassinet, Open

Depth	31"	Depth	31"	Depth	31"	Depth	30 7/8"	Depth	30 7/8"
Width	18 1/4"	Width	18 1/4"	Width	18 1/4"	Width	19 3/4"	Width	19 3/4"
Height	37 3/4"	Height	37 3/4"	Height	37 3/4"	Height	37 1/4"	Height	37 1/4"
Height with Basket	41"	Height with Basket	41"	Height with Basket	41"	Height with Basket	44 1/8"	Height with Basket	44 1/8"