

Pediatric Series Phlebotomy Chairs



Healthcare Furniture

Pediatric Series Phlebotomy Chairs



Bring colour, comfort, and confidence to pediatric blood draw environments with phlebotomy chairs designed to appeal to children while providing the durability and functionality clinicians rely on. Vibrant graphics are embedded directly into the easy-clean laminate to help maintain their bright appearance through routine disinfecting, while the flip arm, right armrest, thickly padded seat and backrest, and medical-grade vinyl upholstery provide a welcoming and supportive experience for both young patients and adults.



66011

Pediatric Blood Draw Chair

Depth	27"
Width	31"
Height	42"
Seat Depth	18"
Seat Width	17.5"
Seat Height	26"
Arm Height	33-39"
Weight Capacity	400 lbs



66099

Extra-Wide Pediatric Blood Draw Chair

Depth	27"
Width	46"
Height	42"
Seat Depth	18"
Seat Width	35"
Seat Height	26"
Arm Height	33-39"
Weight Capacity	400 lbs

Pediatric Series Phlebotomy Chairs



OPTIONS



Aquarium Arm



Arctic Circle Arm



Spring Garden Arm

