

Power Series Phlebotomy Chairs



Healthcare Furniture

Power Series Phlebotomy Chairs



The Power Series phlebotomy chairs combine a comfortable padded seat with convenient electric height adjustment, making it easy for clinicians to position patients at a suitable working height.



6341

Powered Phlebotomy Chair

Depth	30"
Width	32"
Seat Depth	18"
Seat Width	21"
Seat Height	21.5-29.5"
Weight Capacity	375 lbs



6342

Powered Phlebotomy Chair with Drawer

Depth	30"
Width	38.75"
Seat Depth	18"
Seat Width	21"
Seat Height	21.5-29.5"
Weight Capacity	375 lbs



6351

Extra Wide Powered Phlebotomy Chair

Depth	30"
Width	40"
Seat Depth	18"
Seat Width	30"
Seat Height	21.5-29.5"
Weight Capacity	375 lbs

Power Series Phlebotomy Chairs



6352

Extra Wide Powered Phlebotomy Chair with
Drawer

Depth	30"
Width	48"
Seat Depth	18"
Seat Width	30"
Seat Height	21.5-29.5"
Weight Capacity	375 lbs



6361

Bariatric Powered Phlebotomy Chair

Depth	32"
Width	45"
Seat Depth	20"
Seat Width	35"
Seat Height	21.5-29.5"
Weight Capacity	850 lbs



6362

Bariatric Powered Phlebotomy Chair with
Drawer

Depth	32"
Width	52"
Seat Depth	20"
Seat Width	35"
Seat Height	21.5-29.5"
Weight Capacity	850 lbs

