

Rise Bassinet

 **statmedical**
we take care



Durable, Clean, and Safe

The Rise height-adjustable bassinet showcases seamless surfaces that are easy to clean and built from strong KYDEX® materials that offer proven durability. Its adjustable height—thanks to the patented height adjustment columns—makes it a favorite for Users of all heights and moms recovering from C-sections. Its clean, stainless steel push handle and removable drawer are essential to infection prevention.



Features

- Patented aluminum columns with added strength & stability
- Pneumatic spring with a 5-year warranty
- Impact-resistant PETG plastic tub
- KYDEX® thermofilm bassinet and pass-thru drawer
- Four height adjustment handles that can be used independently
- Optional drop leaf surfaces provide extra space
- Trendelenburg and Reverse Trendelenburg positioning
- Infant mattress included

Options

- Two drawers
- Optional drop-leaf work surfaces for either or both sides
- Slide-out work surface
- Available in six colours

	Depth	Surface Depth	Width	Surface Width	Surface Height
M06	22"	21.25"	37.5"	36"	30.5" - 40.25"

