

Rise

Height Adjustable Bassinet



Healthcare Furniture

Rise

Height Adjustable Bassinet



The Rise Height-Adjustable Bassinet is designed to support Canadian nurses with durable, easy-to-clean KYDEX surfaces and flexible height adjustment. Patented height-adjustment columns and accessible handles accommodate different users and mothers recovering from C-sections, while the base panel allows closer access to patients. A stainless steel push handle and removable drawer tub also support effective infection prevention.



M06-0022

Height Adjustable Bassinet

Depth	22"
Width	35"
Height	30.5-40.25"
Surface Depth	19.5"
Surface Width	34.75"
Drawer Dimensions	22.5 x 13.25 x 6"



M06-0020

Height Adjustable Bassinet, Single Drop Leaf

Depth	22"
Width	35"
Height	30.5-40.25"
Extended Width	44.75"
Surface Depth	19.5"
Surface Width	34.75"
Drawer Dimensions	22.5 x 13.25 x 6"



M06-0021

Height Adjustable Bassinet, Dual Drop Leaf

Depth	22"
Width	35"
Height	30.5-40.25"
Extended Width	54.75"
Surface Depth	19.5"
Surface Width	34.75"
Drawer Dimensions	22.5 x 13.25 x 6"

Rise

Height Adjustable Bassinet

 **statmedical**
we take care

Kydex



American
Natural



Aria



Blanca Elme



Bright White



Eiderdown



Summer Flame



Walnut

