

Rise II

Height Adjustable Bassinet



Healthcare Furniture

Rise II

Height Adjustable Bassinet



The Rise II Bassinet features convenient height adjustment with levers located at all four corners for easy access and operation. Its modern design brings a warm, home-like feel to the care environment while also providing ergonomic support to help caregivers work more comfortably.



M09-0042

Single Drawer

Depth	22"
Width	35"
Height	30.5-40.25"
Surface Depth	21.25"
Surface Width	36"
Drawer Dimensions	22.5 x 13.25 x 6"



M09-0022

Two Drawers

Depth	22"
Width	35"
Height	30.5-40.25"
Surface Depth	21.25"
Surface Width	36"
Drawer Dimensions	22.5 x 13.25 x 6"

Rise II

Height Adjustable Bassinet

Kydex



American
Natural



Aria



Blanca Elme



Bright White



Eiderdown



Summer Flame



Walnut

