

Roam

Adjustable Workstation



Built to Last and Fit Your Height

The Roam Mobile Height Adjustable Workstation is a durable, easy-to-clean surface that brings comfort and convenience to any room. Its powerful pneumatic column offers smooth height adjustment without sacrificing surface strength, making it an excellent choice for caregivers wherever they need it.



Features

- Patented Aluminum Height adjustment column that offers superior strength, stability, and smooth height adjustment
- Pneumatic Spring with a 5-Year warranty
- Top utilizes KYDEX thermofilm, 4x thicker than standard thermofilm

Options

- Solid colour
- Wood grain finish
- Drawer
- Locking Drawer

Depth

H00

20"

Width

23.5"

Height

29" - 46"

