

SC Series

Phlebotomy Chair



Healthcare Furniture

SC Series

Phlebotomy Chair



The SC Series of phlebotomy chairs offers a practical and affordable seating solution for a range of budgets. With a minimal footprint, comfortable padded seating, and padded arms available in two sizes, the series provides flexible functionality for a variety of healthcare settings.



6050-U

Standard, Padded Arms



6059-U

Standard, Drawer Arm



6006-U

Extra-Wide, Padded Arms



6069-U

Extra-Wide, Drawer Arm

Depth	25"	Depth	26.5"	Depth	25"	Depth	26.5"
Width	28"	Width	36"	Width	43"	Width	51.5"
Height	18"	Height	18"	Height	18"	Height	18"
Seat Depth	18"	Seat Depth	18"	Seat Depth	18"	Seat Depth	18"
Seat Width	17.5"	Seat Width	17.5"	Seat Width	30"	Seat Width	30"
Seat Height	18"	Seat Height	18"	Seat Height	18"	Seat Height	18"
Arm Height	27-33"	Arm Height	27-33"	Arm Height	27-33"	Arm Height	27-33"
Weight Capacity	300 lbs	Weight Capacity	300 lbs	Weight Capacity	300 lbs	Weight Capacity	300 lbs