

SC500

Sleeper Recliner



Simple and Comfortable

The classic, straightforward design of the SC500 Sleeper Recliner offers infinite reclining possibilities. The leg rest is operated by a simple lever, and the back by a finger pull that lets the chair recline all the way back to a fully flat position. The casters ensure safe transportation, and they lock in place for further security.



SC500

Sleeper Recliner



Features

- Composite foam seat with dynamic suspension
- Shielded mechanism for safety
- Infinite adjustable back
- Lay-flat capability
- 500 lbs capacity
- Heavy-duty mechanism covered by 12-year warranty

Options

- Trendelenburg, motorized
- Continuous footrest
- Glider
- Wood or polyurethane arm caps
- Central locking or directional locking casters
- 3" or 5" casters
- Heated seat and lumbar
- Straight or flared back
- Folding side or profile tablet
- Ergotron arm
- Drainage bag holder

SC500-11

Depth 38"
Width 27"
Seat Width 24"
Height 44"
Seat Height 18"
Weight Capacity 350 lbs



Three-position
Sleeper Recliner

SC500-51M

Depth 37"
Width 32"
Seat Width 24"
Height 45"
Seat Height 18"
Weight Capacity 500 lbs



Multilock Sleeper
Recliner

SC500-51

Depth 38"
Width 32"
Seat Width 23"
Height 45"
Seat Height 18"
Weight Capacity 500 lbs



Multilock Sleeper
Recliner

SC500-52

Depth 37"
Width 32"
Seat Width 24"
Height 45"
Seat Height 18"
Weight Capacity 500 lbs



Multilock Sleeper
Recliner with Trendelenburg

SC500-53

Depth 37"
Width 32"
Seat Width 24"
Height 46"
Seat Height 18"
Weight Capacity 500 lbs



Gliding Multilock
Sleeper Recliner



statmedical
we take care