

SC500

Sleeper Recliner



Healthcare Furniture

SC500

Sleeper Recliner



The SC500 Sleeper Recliner features an independently operated footrest controlled by a lever and an infinite-position backrest operated by a finger pull. The backrest can recline fully flat, providing added flexibility for rest and comfort.



SC500-11



SC500-51



SC500-51M



SC500-52



SC500-53

Three Position Sleeper		Multilock Sleeper Recliner		Motorized Sleeper Recliner		Multilock Sleeper Recliner with Trendelenburg		Gliding Sleeper Recliner	
Depth	39"	Depth	38"	Depth	38.5"	Depth	37"	Depth	37"
Width	30"	Width	32"	Width	33"	Width	31"	Width	32"
Height	44"	Height	45"	Height	45"	Height	46"	Height	45"
Seat Depth	22.5"	Seat Depth	24"	Seat Depth	23"	Seat Depth	23"	Seat Depth	23"
Seat Width	24"	Seat Width	23"	Seat Width	24"	Seat Width	24"	Seat Width	24"
Seat Height	19"	Seat Height	18"	Seat Height	20"	Seat Height	18"	Seat Height	18"
Reclined Depth	68"	Reclined Depth	72"	Reclined Depth	73"	Reclined Depth	72"	Reclined Depth	72"
Arm Height	25.75"	Arm Height	26"	Arm Height	26"	Arm Height	27"	Arm Height	26"
Weight Capacity	350 lbs	Weight Capacity	550 lbs	Weight Capacity	550 lbs	Weight Capacity	550 lbs	Weight Capacity	550 lbs

SC500

Sleeper Recliner

OPTIONS



Moisture Barrier



Drainage Bag Hook



Rolling Foot Tray



IV Pole



Footrest Guard



Fall Prevention Footrest



Push Bar



Adjustable Head Bolster



Flared Back

SC500

Sleeper Recliner

Polyester



Wood

