

SC910

Sleep Chair



Safety-Focused Comfort

An infection control favorite, the SC-910 Sleep chair boasts a simple design with increased comfort. Its metal frame, protected from spills and infestations, is excellent for bariatric patients, as it can hold strong at a 750lb capacity. This chair's design is intuitive and simple, as it has no mechanical parts but can still transition from chair to bed easily at the pull of a lever. It is also easy to handle, maintain, and safely transport.



Features

- A minimalist mechanism comprised of one stainless steel, highly durable spring
- Durable one-piece frame construction
- Easy to open and close
- Locking casters
- Wall-saver design
- Replaceable parts

Options

- Inside Sealed Seams that protect the components and resist infestations
- Contrasting Upholstery for a more personalized appearance

	Depth	Width	Seat Width	Height	Seat Height	Weight Capacity
SC910	36"	35"	27"	38.5"	18"	750 lbs

