

SMART

Rollboards



Minor Equipment

SMART Rollboards



The SMART Rollboard allows for smooth, no-lift patient transfers between surfaces, reducing physical strain on caregivers while providing a safe and comfortable experience for patients.



4100

Standard Foldable Rollboard

Depth	70"
Width	20"
Weight	2.7 kg



4500

MRI/CT Foldable Rollboard

Depth	70"
Width	16"
Weight	2.5 kg



4460

Surgiboard Wide Rollboard

Depth	48"
Width	20"
Weight	1.5 kg



4400

Surgiboard Long Narrow Rollboard

Depth	48"
Width	16"
Weight	1.5 kg

SMART Rollboards



4450

Surgiboard Short Narrow Rollboard

Depth	40"
Width	16"
Weight	1.3 kg



4300

Mini Rollboard

Depth	40"
Width	20"
Weight	1.5 kg



4350

Mini-Gyn Rollboard

Depth	36"
Width	20"
Weight	1.4 kg

