

# SO-400

Station One Bassinet



# Soothing Safety, Smooth Transport

The Station One Bassinet line is designed to inspire ease of mind in patients and their parents alike. The SO-400 is available in healthcare-grade laminate, thermofoil, kydex, antimicrobial, or stainless steel construction. Its self-closing drawers prevent them from opening while moving. The structure is supported by 5" casters designed for smooth, safe transport.



## Features

- Soft close drawer slides prevent accidental opening
- Stainless steel handles
- 5" casters, two locking

## Options

- Laminate, thermofoil, kydex, antimicrobial, or stainless steel construction
- Slide-out writing surface
- Drawer dividers
- Chart holder
- Electronic and key lock options
- Two standard configurations, custom available

Depth

Width

Height

SO-412

21"

34"

37.5"

SO-420

21"

34"

37.5"



SO-412



SO-420