

SO3014-ST

Hydraulic Lift Stool



Healthcare Furniture

SO3014-ST

Hydraulic Lift Stool



The SO3014-ST Hydraulic Stool features a 3.5-inch-thick high-resilience foam seat designed to maintain its original shape while providing long-term comfort. The stool adjusts from 24 to 31 inches in height, and the rotating seat can be locked securely in position by raising the hydraulic control bar.

Depth	23"
Width	23"
Height	36-43"
Seat Width	16"
Seat Height	24-31"
Weight Capacity	400 lbs