

SSC-211

Hydraulic Stool



Healthcare Furniture

SSC-211

Hydraulic Stool



The SSC-211 Hydraulic Surgeon Chair is ideal for surgery, anesthesiology, ENT, and reconstructive procedures. Engineered using the STERIL-GARD philosophy, its design helps minimize or eliminate areas where dirt, fluid, and bacterial debris can become trapped. This makes the chair easier to clean, supports longer-lasting performance, and helps reduce the transmission of bacteria.

Seat Diameter	16"
Base Diameter	21"
Height	22.75-30.25"
Weight Capacity	300 lbs

SSC-211

Hydraulic Stool

OPTIONS

