

SSC-212

Hydraulic Stool



Healthcare Furniture

SSC-212

Hydraulic Stool



The SSC-212 Hydraulic Surgeon Chair is designed for surgery, anesthesiology, ENT, and reconstructive procedures. Built using the STERIL-GARD philosophy, its construction helps minimize or eliminate areas where dirt, fluid, and bacterial debris can become trapped. This makes the chair easier to clean, supports longer-lasting performance, and helps reduce the transmission of bacteria.

Base Diameter	21"
Height	20.75-28.25"
Weight Capacity	300 lbs

SSC-212

Hydraulic Stool

OPTIONS



Moisture Barrier



Stationary Glides



Adjustable Arms

