

SSC-220

Chrome Two Step Stool



Minor Equipment

SSC-220

Chrome Two Step Stool



The SSC-220 Chrome Two Step Stool features a durable chrome finish, extra-wide and deep steps, and a reinforced base for added strength and stability. Its first step measures 10" high x 12.25" x 17.25", with a 350 lb weight capacity under normal use. The handrail can be mounted on either side for convenient access.



SSC-220

Two-Step Chrome Step Stool	
Depth	27"
Width	17.25"
Platform Height	10/17.25"



SSC-220H

Two-Step Chrome Step Stool with Handrail	
Depth	27"
Width	17.25"
Platform Height	10/17.25"
Handle Height	43"