

Solace

Phlebotomy Chair



Healthcare Furniture

Solace

Phlebotomy Chair



The Solace features two infinitely adjustable armrests, an upholstered seat and back, and non-marring levelling feet for added stability. A flip-up contoured urethane armrest provides comfort, easy access, and added patient security. The chair is available with or without a storage cabinet.



2572

Phlebotomy Chair with Storage Drawer

Depth	27"
Width	35.25"
Height	44.5"
Seat Depth	18.5"
Seat Width	18.75"
Seat Height	22.5"
Arm Height	27.75-33.75"
Weight Capacity	300 lbs



2573

Standard Phlebotomy Chair

Depth	27"
Width	29"
Height	44.5"
Seat Depth	18.5"
Seat Width	18.75"
Seat Height	22.5"
Arm Height	27.75-33.75"
Weight Capacity	300 lbs