

Spirit

Phlebotomy Chair



Healthcare Furniture

Spirit

Phlebotomy Chair



The Spirit features two infinitely adjustable armrests, a durable resin seat, and non-marring levelling feet for added stability. Its flip-up contoured urethane armrest provides comfort, easy access, and added safety for the user. The chair is available with or without a storage cabinet.



2570

Phlebotomy Chair with Storage Drawer

Depth	20.5"
Width	35.25"
Height	36"
Seat Depth	15"
Seat Width	18"
Seat Height	20.5"
Arm Height	26.5-32.5"
Weight Capacity	300 lbs



2571

Standard Phlebotomy Chair

Depth	20.5"
Width	29"
Height	36"
Seat Depth	15"
Seat Width	18"
Seat Height	20.5"
Arm Height	26.5-32.5"
Weight Capacity	300 lbs