

Split-Top Overbed Table



Healthcare Furniture

Split-Top Overbed Table



The Split-Top Overbed Table provides added functionality for both patients and caregivers. Its secondary top slides out from either side of the upper surface, while the lower top is released with a tab and can be moved easily with one hand. When the tab is pressed, the two tops move in opposite directions to help balance the load.



ST170

C Base Overbed Table

Depth	18"
Width	32"
Height	29.5-44.5"
Base Height	3.25"
Weight	80 lbs
Weight Capacity	25 lbs



ST170H

H Base Overbed Table

Depth	18"
Width	32"
Height	29.5-44.5"
Base Height	3.25"
Weight	80 lbs
Weight Capacity	25 lbs

Split-Top Overbed Table

OPTIONS



Laminate




Aria



Candlelight



Chestnut



Chocolate



Dolce Vita



Grenada



Hardrock Maple



Havana



Palomino



Pebble



Sand



White Chocolate