

Station One

Double Basin Solution Stand



Minor Equipment

Station One

Double Basin Solution Stand



The Station One Double Basin Solution Stand features durable, fully welded stainless steel construction and rubber-covered ball-bearing casters for smooth mobility. Available with either a lower shelf or support bracing, it includes two 8.5-quart bowls for convenient use in healthcare environments.



H-Brace

MODEL	DEPTH	WIDTH	HEIGHT
SO1002	15"	33"	34"



Shelf

MODEL	DEPTH	WIDTH	HEIGHT
SO1003	15"	33"	34"