


Station One


Free Access Arm and Hand Table with 2" Pad



Minor Equipment

Station One

Free Access Arm and Hand Table with 2" Pad



The SO301 Station One Surgical Table provides unobstructed access to the surgical site with no leg interference, allowing full C-arm and X-ray positioning. Its radiolucent centre section supports clear imaging, while the 12" x 30" operating surface securely attaches to standard OR tables using the Posi-Lock system. Ideal for fistula creation, angiograms, and angioplasty procedures, it includes a 2" comfort pad.

Depth	30"
Width	12"