

# Station One

Hourglass Arm and Hand Table with Leg



## Minor Equipment

# Station One

## Hourglass Arm and Hand Table with Leg



The SO320 Station One Hourglass Arm and Hand Table with Leg features a radiolucent, hourglass-shaped surface for excellent C-arm imaging and surgeon comfort. Its 15" x 32" operating surface expands to 24" at the ends, with reinforced edges and integrated Posi-Lock clamps for secure attachment to OR tables. Available with single or double foot configurations, it offers up to 12" of leg adjustment, an optional Add-A-Rail for accessories, and an included 2" pad.



### SO320

Single Leg

|       |           |
|-------|-----------|
| Depth | 32"       |
| Width | 24/15/24" |



### SO321

Double Leg

|       |           |
|-------|-----------|
| Depth | 32"       |
| Width | 24/15/24" |