

Station One

Patient Slider Boards



Minor Equipment

Station One

Patient Slider Boards



The Station One Slider Board is constructed from durable, high-density polyethylene with smooth surfaces on both sides for easy patient transfers and cleaning. Integrated handholds around the board provide secure grip options, helping caregivers perform safer, more efficient patient transfers.



SO122-W

Bariatric Full-Length Slider Board

Depth	72"
Width	32"



SO122

Full-Length Slider Board

Depth	72"
Width	22"



SO121

Half-Size Slider Board

Depth	33"
Width	22"